



Brand: Honeywell Riedel-de Haën
Product: 24103

Ethanol puriss., absolute, sales not in Germany, 99.8% (GC)

Parameter	Specification	Units
assay (GC)	min. 99.8	Vol.%
density (D 20/20)	0.790-0.791	
non-volatile matter	max. 0.01	%
water (Karl Fischer)	max. 0.2	%
free acid (as CH ₃ COOH)	max. 0.001	%
aldehydes (as CH ₃ CHO)	max. 0.005	%
free alkali (as NH ₃)	max. 0.0005	%
methanol (GC)	max. 0.1	%
S-compounds (as S)	max. 0.2	ppm