



AS 4000.2 Fixing clip attachment

For shaking of flasks, Erlenmeyer flasks and bottles with a round crosssection (without fixing clips).

Dimensions (W x D x H): 470 x 440 x 25 mm

Number of fixing clips (volume): 50 x AS 2.1 (25 ml), 48 x AS 2.2 (50 ml), 25 x AS 2.3 (100 ml), 16 x AS 2.4 (250 ml), 12 x AS 2.5 (500 ml)

Set-up plate: 430 x 430 mm

Weight: 2.650 g

Accessories: AS 2.1 Fixing clips, AS 2.2 Fixing clips, AS 2.3 Fixing clips, AS 2.4 Fixing clips, AS 2.5 Fixing clips