


Fiche technique

Mélange / Agitateurs-secoueurs



Trayster basic

Overhead shaker with vertical rotation for smooth yet effective mixing at a fixed speed of 20 rpm. Ideal for biological samples, particularly blood. Suitable for separating and mixing processes of powdered and liquid samples in Eppendorf and Greiner tubes up to 50 ml.

Easily removable attachments (three attachments fit on one drive) provide for quick cleaning in case of sample spills.

- 1) Different attachments allow a variety of applications.
- 2) Usable at 4 to 50 °C in incubators.
- 3) Suitable for continuous operation.

Accessoires: TR 1 Tube insert, TR 2 Tube insert, TR 3 Tube insert, TR 4 Tube insert

Données techniques

Type de mouvement	hélice
Capacité admissible (support inclus) [kg]	2
Durée de fonctionnement admissible [%]	100
Vitesse min. [rpm]	0
Vitesse fixe [rpm]	20
Affichage de la vitesse	aucun
Minuterie	non
Affichage de la minuterie	aucun
Type de fonctionnement	mode continu
Fonction " touch"	non
Unité superposable	non
Dimensions (L x P x H) [mm]	460 x 180 x 140
Poids [kg]	1.6
Plage de température du milieu admise [°C]	4 - 50
Humidité relative admissible [%]	80
Protection selon DIN EN 60529	IP 21
Interface numérique USB	non
Interface numérique RS 232	non
Sortie analogique	non
Tension [V]	100 - 240
Fréquence [Hz]	50/60
Tension continue [V=]	24
Consommation électrique max. [mA]	1000
Ident. N°.	0004005000