



AS 501.4 Fixing clip support

For shaking flasks, Erlenmeyer flasks and pear-shaped flasks (without fixing clips).

Dimensions (W x D x H): 475 x 460 x 95 mm

Capacity: 110 x AS 2.1 (25 ml) fixing clips (volume): 55 x AS 2.2 (50 ml), 35 x AS 2.3 (100 ml), 16 x AS 2.4 (250 ml), 12 x AS 2.5 (500 ml), 8 x AS 2.6 (1000 ml)

Weight: 2.640 g