



PT 1000.60 Temperature sensor, stainless steel

Made of stainless steel, for use with RCT basic, RET basic and C-MAG HS digital series.

Depth of immersion: 230 mm

Diameter: 3 mm

Material: AISI 316 Ti

Technical Data

Probe type	Immersion probe
Sensor type	PT1000 (DIN IEC 751 Class A)
Accuracy of temperature measurement [K]	± (0,15 + 0,2% o.Mv.)
Response time $t_{9/10}$ [sec]	10
Temperature measuring range min. [°C]	-50
Temperature measuring range max. [°C]	+400
Material of probe tube	stainless steel 1.4571
Diameter of probe tube [mm]	3
Length of probe tube [mm]	230
Immersion depth min. [mm]	25
Cable length [mm]	1000
Probe connection	DIN 45322
Ident. No.	0003516800