

Technical Datasheet

Article No. 2513

1-Butanol p. A.
min. 99.5 % CH₃(CH₂)₃OH

For laboratory use.

Parameter	Wert
Appearance	colourless liquid
Boiling point	117 °C
Density	0.81 g/ml
Mol.-Weight	74.12 g/mol
Non volatile substances	max. 10 mg/kg
Free acid (as CH ₃ COOH)	max. 20mg/kg
Refractive index (20 °C)	1.397–1.401
Water (KF)	max. 1000 mg/kg
Colour (Hazen)	max. 10