

Azlon® Polycarbonate bottle, Boston round

Classic, round design

Featuring a classic, round design, the Azlon® boston round polycarbonate bottle is great for the preparation and storage of biological media and is autoclavable at 121°C.

- Highly transparent polycarbonate media bottles with linerless polypropylene cap
- Suitable for repeated sterilisation by autoclaving at 121°C
- Ideal for the preparation and storage of biological media, reagents and aqueous solutions
- Leakproof
- Classic round design
- Heavy duty and durable in use



Product Specification

Capacity ml	Diameter mm	Height mm	Neck ID mm	Pack quantity
1000	97	216	28	12