

**ClearLine serology pipette**

- Crystal polystyrene
- Sterile (gamma irradiation)
- DNase and RNase free
- Non-pyrogenic
- Color coded for better identification
- Color: Blue
- Cotton pipette: yes
- Capacity: 5 ml
- Total length: 348.15 mm
- Graduation: Bidirectional, 0.1 ml
- Over-graduation: 2.5 ml
- Accuracy: +/- 2%
- Quantity per bag: 50
- Quantity per case: 200

	Lenght	Internal diameter	External diameter
Tip	37,36 mm	1,57 mm	3,35 mm
Tip	N/A	8,03 mm	6,01 mm
Body	311,15 mm	6,01 mm	8,03 mm

**Units : 200**