

Clearline Micropipette

Features :

- ✓ Volume locking prevents accidental counter shift
 - ✓ Ergonomic : contoured handle to adapt to the hand, finger support to reduce wrist fatigue
 - ✓ Ejector with infinitely variable length adjustment makes it easier for different tips
 - ✓ Low pipette weight
 - ✓ Soft springs system provides comfort at work
 - ✓ Fully autoclavable
 - ✓ UV resistant
 - ✓ Calibration key included
 - ✓ Each instrument passes stringent validation and control procedures in accordance with EN ISO 8655
-
- ✓ Volume : 10 - 100 μ l
 - ✓ Color : yellow
 - ✓ Channel : singlechannel



Units : 1 micropipette