

## Clearline Micropipette

### Features :

- ✓ Volume locking prevents accidental counter shift
  - ✓ Ergonomic : contoured handle to adapt to the hand, finger support to reduce wrist fatigue
  - ✓ Ejector with infinitely variable length adjustment makes it easier for different tips
  - ✓ Low pipette weight
  - ✓ Soft springs system provides comfort at work
  - ✓ Fully autoclavable
  - ✓ UV resistant
  - ✓ Calibration key included
  - ✓ Each instrument passes stringent validation and control procedures in accordance with EN ISO 8655
- 
- ✓ Volume : 5 - 50  $\mu$ l
  - ✓ Color : yellow
  - ✓ Channel : 8 channels



**Units : 1 micropipette**