


CleaTime II timer

Specification :

- ✓ Up to 2 counters and / or 2 decounters
- ✓ Clock, alarm, snooze functions (alarm delayed by 5min)
- ✓ Fast increment
- ✓ Memory of the time programmed by channels
- ✓ 2 alarms 1mm
- ✓ Capacity : 24h per second
- ✓ Time exceeded
- ✓ Screen : 40 x 17 mm
- ✓ Number : 15 mm
- ✓ Dimension : 67 x 55 x 20 mm
- ✓ With crutch, clamp and magnet
- ✓ Power supply: A76 / LR44 cell

**Unités :**

1 timer