

SLS Lab Pro Unstirred Water Baths



The SLS Lab Pro Range of Unstirred Water Baths 2-28L provides class leading performance and reliability for scientific, laboratory and industrial applications for consistent temperature control.

Features

PID temperature control, increments 0.1°C

Bold LED temperature display

Stainless steel tank with concealed heaters and sensors

Painted surfaces feature antibacterial finish

Safety features, thermal cut out and over temperature alarm protecting lab work/samples

Factory temperature calibration, fully tested

Accessories available, enquire for details - Lid included

Construction : Stainless steel 304 spec tank. Painted surfaces feature “anti-bacterial technology”

Temperature range °C : Ambient 5 to 99, PID digital LED display showing actual or set temperature, resolution 0.1°C - Stability : $\pm 0,2$ °C

Safety : Over temperature thermal cut out

Voltage : 230 V – Energy use : 0,2 kWh @ 37°C

Warranty : 3 years

Overall Dims (W x D x h) : 332 x 185 x 135 mm – Tank Dims (W x D x h) : 300 x 150 x 65 mm

Tank max capacity : 2 L

Drain : no

Weight : 5,5 kg

Unit : 1 Water Bath