

## SLS Lab Pro Unstirred Water Baths



The SLS Lab Pro Range of Unstirred Water Baths 2-28L provides class leading performance and reliability for scientific, laboratory and industrial applications for consistent temperature control.

### **Features**

PID temperature control, increments 0.1°C

Bold LED temperature display

Stainless steel tank with concealed heaters and sensors

Painted surfaces feature antibacterial finish

Safety features, thermal cut out and over temperature alarm protecting lab work/samples

Factory temperature calibration, fully tested

Accessories available, enquire for details - Lid included

Construction : Stainless steel 304 spec tank. Painted surfaces feature “anti-bacterial technology”

Temperature range °C : Ambient 5 to 99, PID digital LED display showing actual or set temperature, resolution 0.1°C - Stability :  $\pm 0,2$  °C

Safety : Over temperature thermal cut out

Voltage : 230 V – Energy use : 0,2 kWh @ 37°C

Warranty : 3 years

Overall Dims (W x D x h) : 361 x 332 x 290 mm – Tank Dims (W x D x h) : 325 x 300 x 150 mm

Tank max capacity : 14 L

Drain : no

Weight : 8 kg

### **Unit : 1 Water Bath**