

## SLS Lab Pro Unstirred Water Baths



The SLS Lab Pro Range of Unstirred Water Baths 2-28L provides class leading performance and reliability for scientific, laboratory and industrial applications for consistent temperature control.

### **Features**

- PID temperature control, increments 0.1°C
- Bold LED temperature display
- Stainless steel tank with concealed heaters and sensors
- Painted surfaces feature antibacterial finish
- Safety features, thermal cut out and over temperature alarm protecting lab work/samples
- Factory temperature calibration, fully tested
- Accessories available, enquire for details - Lid included

Construction : Stainless steel 304 spec tank. Painted surfaces feature “anti-bacterial technology”  
Temperature range °C : Ambient 5 to 99, PID digital LED display showing actual or set temperature, resolution 0.1°C - Stability :  $\pm 0,2$  °C  
Safety : Over temperature thermal cut out  
Voltage : 230 V – Energy use : 0,55 kWh @ 37°C  
Warranty : 3 years  
Overall Dims (W x D x h) : 538 x 332 x 290 mm – Tank Dims (W x D x h) : 500 x 300 x 150 mm  
Tank max capacity : 22 L  
Drain : yes  
Weight : 11 kg

### **Unit : 1 Water Bath**