

ClearLine syringe filters



The CA syringe filters have a hydrophilic membrane that allows the filtration of aqueous solutions, especially in biology. The sterile syringe filters are suitable for cell culture applications, for proteins or enzymes samples

Features :

- Hydrophilic membrane
- Low protein adsorption
- Limited solvent resistance

Specifications :

- Material : CA - Cellulose acetate
- Color code : Orange
- Pore size : 0,45 µm
- Diameter : 25 mm
- Water flow (mL/mn@10psi) : 92-116
- Bubble point (psi) : 36-39
- pH : 3-7
- Sterile : Yes, gamma-ray methods
- Prefilter : -

Units :

50 syringe filters