

### ClearLine syringe filters



Naturally hydrophilic, with a good mechanical strength and strong absorption, the NYL syringe filters are suitable for the filtration of HPLC samples, filtration and clarification of solvents

**Features :**

- Hydrophilic membrane
- Good chemical resistance

**Specifications :**

- Material : NYL - Polyamide
- Color code : Yellow
- Pore size : 0,45 µm
- Diameter : 13 mm
- Water flow (mL/mn@10psi) : 47-55
- Bubble point (psi) : 27-33
- pH : 3-14
- Sterile : Yes, gamma-ray methods
- Prefilter : -

**Units :**

50 syringe filters