

ClearLine syringe filters



PVDF syringe filters are hydrophobic and offer an excellent chemical resistance against most organic solvents and aggressive liquids. They are suitable for air and gas filtration.

Features :

- Hydrophobic membrane for protein-laden solutions
- Very low protein adsorption
- Good solvent resistance

Specifications :

- Material : PVDF - Polyvinylidene
- Color code : Beige
- Pore size : 0,45 µm
- Diameter : 33 mm
- Water flow (mL/mn@10psi) : 70-105
- Bubble point (psi) : 11-23
- pH : 1-14
- Sterile : No
- Prefilter : -

Units :

100 syringe filters