

### ClearLine syringe filters



RC syringe filters are recommended for the direct filtration of aqueous solutions or of any type of organic solvent. Their benefit is to offer a great capacity of particles extraction.

**Features :**

- Hydrophilic membrane
- Broad chemical compatibility
- Good particle extraction capacity

**Specifications :**

- Material : RC - Regenerated cellulose
- Color code : Blue
- Pore size : 0,45  $\mu\text{m}$
- Diameter : 25 mm
- Water flow (mL/mn@10psi) : 53-74
- Bubble point (psi) : 45-49
- pH : 3-12
- Sterile : Yes, gamma-ray methods
- Prefilter : -

**Units :**

50 syringe filters