

### ClearLine syringe filters



The CA syringe filters have a hydrophilic membrane that allows the filtration of aqueous solutions, especially in biology. The sterile syringe filters are suitable for cell culture applications, for proteins or enzymes samples

**Features :**

- Hydrophilic membrane
- Low protein adsorption
- Limited solvent resistance

**Specifications :**

- Material : CA - Cellulose acetate
- Color code : Orange
- Pore size : 0,45  $\mu\text{m}$
- Diameter : 25 mm
- Water flow (mL/mn@10psi) : 92-116
- Bubble point (psi) : 36-39
- pH : 3-7
- Sterile : No
- Prefilter : With pre-filter for charged solutions

**Units :**

50 syringe filters