

ClearLine syringe filters



RC syringe filters are recommended for the direct filtration of aqueous solutions or of any type of organic solvent. Their benefit is to offer a great capacity of particles extraction.

Features :

- Hydrophilic membrane
- Broad chemical compatibility
- Good particle extraction capacity

Specifications :

- Material : RC - Regenerated cellulose
- Color code : Blue
- Pore size : 0,45 µm
- Diameter : 25 mm
- Water flow (mL/mn@10psi) : 53-74
- Bubble point (psi) : 45-49
- pH : 3-12
- Sterile : No
- Prefilter : With pre-filter for charged solutions

Units :

50 syringe filters