

ClearLine syringe filters



Naturally hydrophilic, with a good mechanical strength and strong absorption, the NYL syringe filters are suitable for the filtration of HPLC samples, filtration and clarification of solvents

Features :

- Hydrophilic membrane
- Good chemical resistance

Specifications :

- Material : NYL - Polyamide
- Color code : Yellow
- Pore size : 0,22 µm
- Diameter : 25 mm
- Water flow (mL/mn@10psi) : 21-36
- Bubble point (psi) : 46-55
- pH : 3-14
- Sterile : No
- Prefilter : With pre-filter for charged solutions

Units :

50 syringe filters