

ClearLine syringe filters



Hydrophobic and with an excellent chemical stability, the PP syringe filters are suitable for HPLC samples containing weak solids or sterile filtration of samples of small volume.

Features :

- Hydrophobic membrane
- Air and solvent filtration
- High throughput

Specifications :

- Material : PP – Polypropylene
- Color code : Natural
- Pore size : 0,22 µm
- Diameter : 25 mm
- Water flow (mL/mn@10psi) : 252-305
- Bubble point (psi) : -
- pH : 1-14
- Sterile : No
- Prefilter : With pre-filter for charged solutions

Units :

50 syringe filters