

### ClearLine syringe filters



Hydrophobic and with an excellent chemical stability, the PP syringe filters are suitable for HPLC samples containing weak solids or sterile filtration of samples of small volume.

**Features :**

- Hydrophobic membrane
- Air and solvent filtration
- High throughput

**Specifications :**

- Material : PP – Polypropylene
- Color code : Natural
- Pore size : 0,45 µm
- Diameter : 33 mm
- Water flow (mL/mn@10psi) : 550-635
- Bubble point (psi) : -
- pH : 1-14
- Sterile : No
- Prefilter : -

**Units :**

50 syringe filters