


**DATASHEET****PART NUMBER****257297**

## SolidFrame PCR plates

**Description :**

- The rigid semi-skirted frame remains dimensionally stable during thermal cycling, which facilitates reliable robotic handling and ensures that the seal to each well is not compromised during thermal cycling
- Thin-walled tubes made from 100% virgin medical grade polypropylene
- Rigid polycarbonate frame
- Low profile wells have a maximum capacity of 225 µl
- Manufactured to the highest quality standards for your PCR and qPCR applications
- Our 2-component Roche style design delivers improved performance for the lightcycler
- Certified free from Dnase, Rnase, nucleases and human gDNA

**Specification :**

- Number of well : 96
- Frame color : natural
- Well color : white
- Skirt : half-skited
- Style : Roche
- Well volume : 0,1 ml

**Units : 50 plates**