


DATASHEET**PART NUMBER****257301****SolidFrame PCR plates****Description :**

- The rigid semi-skirted frame remains dimensionally stable during thermal cycling, which facilitates reliable robotic handling and ensures that the seal to each well is not compromised during thermal cycling
- Thin-walled tubes made from 100% virgin medical grade polypropylene
- Rigid polycarbonate frame
- Standard profile wells have a maximum capacity of 330 µl
- Manufactured to the highest quality standards for your PCR and qPCR applications
- ABI style 2-component plates are designed for compatibility with most ABI thermal cyclers, from legacy designs to the most recent instruments
- Certified free from Dnase, Rnase, nucleases and human gDNA

Specification :

- Number of well : 96
- Frame color : natural
- Well color : natural
- Skirt : half-skited
- Style : ABI
- Well volume : 0,2 ml

Units : 50 plates