


DATASHEET**PART NUMBER****257306****SolidFrame PCR plates****Description :**

- The rigid skirted frame remains dimensionally stable during thermal cycling, which facilitates reliable robotic handling and ensures that the seal to each well is not compromised during thermal cycling
- Thin-walled tubes made from 100% virgin medical grade polypropylene
- Rigid polycarbonate frame
- Maximum well capacity is 45 µl
- Manufactured to the highest quality standards for your PCR and qPCR applications
- Our fully skirted 2-component PCR plate fits a wide range of instruments and is fully automation compatible
- Certified free from Dnase, Rnase, nucleases and human gDNA

Specification :

- Number of well : 384
- Frame color : black
- Well color : natural
- Skirt : skirted
- Style : -
- Well volume : - ml

Units : 50 plates