


**20µl Clearline Prime filter tips sterile graduated low retention****Features**

- ClearLine® Prime long graduated filter tips 20 µl
- Material: Natural Virgin Polypropylene Tips
- With filter
- Sterile (e-beam)
- Low retention
- Universal compatibility
- Non autoclavable
- Storage: Store in a clean, dry place at room temperature.
- Quality control: Each batch is certified: RNase, DNase-free, free of human DNA, endotoxins and pyrogens
- Traceability: On each rack: a 5-digit batch number, the reference and the volume of the tips as well as their expiry date.

**Conditioning**

- 10 x 96 tips, in hinged rack refill system

**Units :**

960 Tips