


### 1000µl Clearline Prime natural reload tips non sterile



#### Features

- ClearLine® Prime graduated tips 1000 µl
- Material: Natural Virgin Polypropylene Tips
- Without filter
- Non sterile
- Universal compatibility
- Autoclavable
- Storage: Store in a clean, dry place at room temperature.
- Quality control: Each batch is certified: RNase, DNase-free, free of human DNA, endotoxins and pyrogens
- Traceability: On each rack: a 5-digit batch number, the reference and the volume of the tips as well as their expiry date.

#### Conditioning

- 4 x 960 tips

#### Units :

3840 Tips