

**1250 µl ClearLine® Prime natural reload tips non sterile low retention****Features**

- ClearLine® Prime long graduated tips 1250 µl
- Material: Natural Virgin Polypropylene Tips
- Without filter
- Non sterile
- Low retention
- Universal compatibility
- Autoclavable
- Storage: Store in a clean, dry place at room temperature.
- Quality control: Each batch is certified: RNase, DNase-free, free of human DNA, endotoxins and pyrogens
- Traceability: On each rack: a 5-digit batch number, the reference and the volume of the tips as well as their expiry date.

**Conditioning**

- 4 x 960 tips

**Units :**

3840 Tips