

**1250µl Clearline Prime natural reload non sterile****Features**

- ClearLine® graduated 1250 µl unfiltered tips
- Material: Natural Virgin Polypropylene Tips
- Without filter
- Not sterile
- Universal compatibility
- Autoclavable
- Storage: Store in a clean, dry place at room temperature.
- Quality control: Each batch is certified: RNase, DNase-free, free of human DNA, endotoxins and pyrogens
- Traceability: On each rack: a 5-digit batch number, the reference and the volume of the tips as well as their expiry date.

**Conditioning**

- 8 x 96 tips, in hinged rack refill system

**Units :**

768 Tips