


**White PCR plate 96 semi-skirted wells 0.1ml**

- 96-well white PCR plate
- Semi-skirted wells 0.1 ml
- DNase, RNase and gDNA free
- Pyrogen and endotoxin free
- Labeled surface: horizontal numbers (1-12) and vertical letters (A-H)
- Ultra-thin wall technology provides excellent heat transfer effect which leads to maximum amplification
- Flange and tapered tube design for effective sealing to prevent evaporation and cross contamination
- Compatible with most automatic devices
- Dimensions (L x W x H): 128.3mm x 85.5mm x 15.7mm

**Units :**

50 PCR plates