

### ClearLine extraction thimbles

With the appropriate solvents, the extraction thimbles enable to extract specific particles from a solid substance, allowing a more accurate and swifter analysis. They are neutral and fat free and offer a high mechanical strength and an excellent retention capacity. The very high quality of the fibers used for the extraction thimbles allow a great capability of reproduction during the analysis.

Material : Cotton linter fibres

Dimensions : 22 x 80 mm

Thickness : 1 mm

Pore size : 5  $\mu$ m

**Unit :**

25 extraction thimbles