

ClearLine glass microfiber filter with binder

The glass microfiber filters with binder have a lower resistance to temperature (up to 180 °C maximum). The hydrophobic grade [FV27] is suitable for air and gas analysis. At the opposite, the hydrophilic grade [FV29] is adapted to liquid filtration.

Application : Gravimetric analysis

Weight DIN 53104 : 73 g/m²

Thickness : 0,35 mm

Propriété : Hydrophilic

Retention efficiency for 0,3 µm : 0

Retention : 0,6 µm

Dimensions : 125 mm

Units :
50 filters