

### ClearLine filter papers

The standard filter papers are made with 100% cellulose and can have various applications in hospitals/ medical, in research laboratories, in schools/universities, in industry such as chemical, pharmaceutical, cosmetic, water treatment, food...

Application : Filter paper for analytical works in laboratories with medium filtration

Weight DIN 53104 : 73 g/m<sup>2</sup>

Thickness : 0,16 mm

Pore size : 5-13 µm

Filtration speed DIN 53137 : 35 sec

Burst strength : 1,95 kg/cm<sup>2</sup>

Format : Plain discs

Dimensions : 500 mm

**Units :**

**100 filters papers**