

ClearLine filter papers

The standard filter papers are made with 100% cellulose and can have various applications in hospitals/ medical, in research laboratories, in schools/universities, in industry such as chemical, pharmaceutical, cosmetic, water treatment, food...

Application : Filter paper for analytical works in laboratories with medium filtration

Weight DIN 53104 : 73 g/m²

Thickness : 0,16 mm

Pore size : 5-13 µm

Filtration speed DIN 53137 : 35 sec

Burst strength : 1,95 kg/cm²

Format : Folded discs

Dimensions : 100 mm

Units :

100 filters papers