

ClearLine Quantitative analytical filters

The quantitative or ashless filter papers are made with 100% high quality cotton linter fibres, which go through a severe chemical process. They have been washed out with a specific acid treatment and finally cleaned from impurities with demineralized water. This process allows reaching an ash content below 0.01%. These papers are used to count during demanding analysis.

Application : Slow filtration quantitative filter paper

Colour code : Green

Weight DIN 53104 : 74 g/m²

Thickness : 0,16 mm

Pore size : 7-9 µm

Filtration speed DIN 53137 : 100 sec

Format : Plain discs

Dimensions : 70 mm

Unités :
100 filters