

### ClearLine membranes filters

The MCE membranes are made with cellulose acetate and cellulose nitrate fibres, they are naturally hydrophilic, mechanically stable and have a high loading capacity. They are suitable for microbiological analysis, for colonies counting or for pre-filtration and clarification of samples.

Membrane : MCE - Mixed cellulose ester

Pore size (for Ø 47 mm) : 3,00 µm

Water flow @ 10 psi : 210-320 ml/min/cm<sup>2</sup>

Bubble point : 8-11 psi

Diameter : 90 mm

Non sterile

Colour : White, plain

**Units :**

100 membranes