


Tips 200 µl ClearLine® natural bulk**Features**

- Natural Bulk Standard Tips
- Volume 200µl
- Scale 10, 50, 100
- Length 51.44mm
- natural color
- Filter no
- Low Retention no
- Not sterile
- Autoclavable
- Packaging: in bulk
- Material: virgin polypropylene
- Store in a clean, dry place at room temperature
- Each batch is certified: RNase, DNase free, free of PCR inhibitor, human DNA, endotoxins, pyrogens and PCR inhibitor
- Traceability: 5-digit batch number, reference, tip volume and expiration date

Units :

1000 Tips