


Clearline tips 10 μ L Hybridrack ClearLine[®] sterile

- In clear polypropylene
- Store in a clean and dry place, at room temperature
- Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	258361
Packaging	HybridRack: cardboard rack with recyclable tip tray. (10 racks of 96 tips)
Volume	10 μ l
Graduation	2, 5, 10
Length	45,34 mm
Color	Natural
Filter	No
Low Retention	Yes
Sterility	Yes
Autoclavable	Yes



Units :
960 Tips