


**200 µl tip, LTS compatible, graduated, sterile, reload, blue tray**  
**ClearLine**

- Volume : 1000 µl
- Longueur : 98,6 mm
- Packaging : Reload Clearline Plus - 8 x 96
- Compatibility : Compatible with Rainin LTS™ micropipettes
- Tray color : blue
- Low retention : no
- Autoclavable : yes
- Sterile : no
- Filter : no
- Color : natural
- Graduation : 200, 500, 1000
- Material : polypropylene

**Units : 768**