



Shorty serological pipettes

- Manufactured from highly transparent polystyrene (PS) material
- With a filter plug for filtration
- Different color bands
- Clear dual-scale lines for positive and negative readings
- Easy to handle under a cabinet
- Same precision as a standard pipette
- Free from DNA enzymes, RNA enzymes, and pyrogene
- Irradiation sterilization
- Package : Individually packaged in paper/plastic
- Volume : 5 ml
- Graduation : 0,100000000000000001
- Color code : Blue
- Length : 232 mm

**Units :**

200