



# Foam Tipped Applicator - Dry Transport System

Thin flexible foam tip, nylon handle



## Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. This specialty tipped applicator is made with 100ppi medical grade polyurethane foam. It is bonded securely to a tapered acetal handle.

## Dimensions

### Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.00" (152.4 mm)

### Tip

- 0.141" Dia. (3.581 mm)
- 0.547" Length (13.894 mm)

### Handle

- 0.098" Dia. (2.489 mm) tapered to 0.049" (1.321 mm)
- 5.969" Length (152.298 mm)

## Qualities

- Medical grade quality
- Foam has excellent collection and release properties
- Soft, small tip for patient comfort
- Lint and residue free
- Made in USA

## Applications

- Specimen collection
- Diagnostic sampling
- Recovery of microorganisms and fluids

Ref #	Packaging
25-1406 1PF BT Sterile	5/100/1 (5 boxes per case, 100 tubes per box, 1 swab per tube)



**Puritan**<sup>®</sup>