



# FAB-MINI-AP Cotton Tipped Applicator

Mini-tip tip, wood handle



## Description

This DNA controlled absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a polystyrene handle.

The convenient, reclosable, dry transport device is designed for reliable sample and evidence collection. The patented tube design entails a breathable filter label that covers holes to prohibit mold growth during transport. ID labels are included.

## Dimensions

5.095" Overall Length (129.413 mm)

### Tip

- 0.12" Dia. (3.048 mm) tapered to 0.05" (1.27 mm)
- 0.5" Length (12.7 mm)

### Handle

- 0.083" Dia. (2.108 mm)
- 5.032" Length (127.813 mm)

## Qualities

- Medical grade quality
- Produced in a DNA controlled environment
- Reclosable transport tube
- Can be used dry, with IPA or sterile water
- Includes evidence information label
- Tapered mini-tip for reaching into small areas
- Individually wrapped
- Made in USA

## Applications

- DNA testing
- Forensic evidence collection

Ref #	Packaging
25-825 1WCS TT FABUSA Sterile	10/50/1 (10 boxes per case, 50 individually wrapped tubes per box, 1 swab p/tube)



**Puritan**<sup>®</sup>