



# HydraFlock® Flocked Swab

Mini-Tip with Polystyrene Handle - no breakpoint



## Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Patented: US, Canada, Europe, Australia, Japan, China  
Patent pending: Hong Kong

## Dimensions

6.04" Overall Length (153.416 mm)

### Tip

- 0.135" Dia. (3.429 mm)
- 0.800" Length (20.32 mm)

### Handle

- 0.098" Dia. (2.489 mm) tapered to 0.049" (1.245 mm)
- 6" Length (152.654 mm)

## Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

## Applications

- For nasopharyngeal and urethral cell collection
- Rapid diagnostic testing
- Specimen collection
- For collection of bacteria or micro-organisms
- Forensic evidence collection

Ref #	Packaging
3316-H Non-Sterile	10/100 (10 bags of 100 swabs per bag)
25-3316-H Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)

