




# PurFlock Ultra® Flocked Swab

Large Tip with Polystyrene Handle - 80mm breakpoint



## Description

PurFlock's purity make it ideal for PCR, molecular assays, DNA testing, forensic applications as well as direct antigen testing. Its uniquely designed fiber structure accounts for better cell yield in the specimen where sample size is limited.

Patented: US, Canada, Europe, Australia, Japan, China  
Patent pending: Hong Kong

## Dimensions

5.93" Overall Length (150.622 mm)

### Tip

- 0.219" Dia. (5.563 mm)
- 0.68" Length (17.27 mm)

### Handle

- 0.100" Dia. (2.54 mm) tapered to 0.09" (2.286 mm)
- 5.89" Length (149.606 mm)

## Qualities

- Highly absorbent tip
- Soft smooth tip
- Made in USA

## Applications

- For collection of bacteria or micro-organisms
- Specimen collection
- Buccal cell collection
- Forensic evidence collection

Ref #	Packaging
3406-U Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3406-U Sterile	10/50/1 (10 boxes per case, 50 individually wrapped swabs, 1 swab p/pkg)

