




# PurFlock Ultra<sup>®</sup> Flocked Swab

Elongated Tip with Polystyrene Handle - 80 mm breakpoint



## Description

PurFlock's purity make it ideal for PCR, molecular assays, DNA testing, forensic applications as well as direct antigen testing. Its uniquely designed fiber structure accounts for better cell yield in the specimen where sample size is limited.

Patented: US, Canada, Europe, Australia, Japan, China  
Patent pending: Hong Kong

## Dimensions

**5.93" Overall Length (150.622 mm)**

### Tip

- 0.219" Dia. (5.563 mm)
- 0.770" Length (19.558 mm)

### Handle

- 0.1" Dia. (2.54 mm) tapered to 0.09" (2.286 mm)
- 5.89" Length (149.606 mm)

## Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

## Applications

- For Buccal cell collection
- Specimen collection
- Rapid diagnostic testing
- Forensic evidence collection
- DNA testing

| Ref #                 | Packaging  |
|-----------------------|--|
| 3506-U<br>Non-Sterile | 10/100<br>(10 bags per case, 100 per bag)              |
| 25-3506-U<br>Sterile  | 10/50/1<br>(10 boxes of 50 individually wrapped swabs) |

