

Rayon Tipped Applicator

Polystyrene handle



- USP Grade Rayon Tip

- Semi-flexible polystyrene handle

Actual Size



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Description

This specialty swab is made with USP grade rayon fiber that is carded on-site to maintain consistency. This swab is ideal for diagnostic specimen collection, wound treatment and application of medicines. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

Dimensions

6" Overall Length (152.4 mm)

Tip

■ 0.203" Dia. (5.156 mm)

■ 0.688" Length (17.475 mm)

Handle

■ 0.100" Dia. (2.54 mm)

■ 5.875" Length (149.225 mm)

Standards

■ CE marked for medical use

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

■ Certificate of sterility and/or compliance available with each shipment upon request

Qualities

■ Medical grade quality

■ Absorbent tip

■ Soft smooth tip

Applications

Diagnostics / Forensics / Medical

■ Specimen collection swab for testing & screening

■ Cell collection for DNA testing

■ General purpose collection swab

Ref #	Packaging
25-806 1PR <i>Sterile</i>	10/100/1 (10 boxes of 100 individually wrapped packages; 1 per package)
25-806 2PR <i>Sterile</i>	10/100/2 (10 boxes of 100 individually wrapped packages; 2 per package)

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com