




Cotton Tipped OB/GYN Applicator

Paper handle



-USP grade cotton tip

- Paper handle

*Image not actual size

Description

This Puritan OB/GYN swab is made with High Quality USP grade cotton which is carded on-site to maintain consistency. The large, oversized, cotton tip is extra-absorbent.

Dimensions

8.188" Overall Length (207.975 mm)

Tip

- 0.5" Dia. (12.7 mm)
- 1.125" Length (28.575 mm)

Handle

- 0.156" Dia. (3.962 mm)
- 8" Length (203.2 mm)

Qualities

- Extra absorbent, large oversized tip
- Long paper handle
- Medical grade quality
- Made in USA

Applications

- Preparatory to specimen sampling or exams
- General cleaning of large surfaces

Ref #	Packaging
808 COTTON Non-sterile	10/500 (10 boxes per case, 50 per box)



Puritan[®]