

Rayon Tipped Proctoscopic Applicator

Polypropylene handle



- USP grade Rayon tip

- Hollow polypropylene handle

* Not actual size



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Description

Puritan Proctoscopic swabs are made with USP grade rayon fiber that is carded on-site to maintain consistency. The large rayon tip is extra-absorbent.

Dimensions

16.062" Overall Length (407.975 mm)

Tip

- 1.125" Length (28.575 mm)
- 0.5" Diameter (12.7 mm)

Handle

- 15.875" Length (403.225 mm)
- 0.150" Dia. (3.81 mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- FDA GMP's
- Autoclavable to 250°F

Qualities

- Extra-absorbent large tip
- Multi use swab
- Medical grade quality

Applications

Medical / Diagnostics / Veterinary

- OB/GYN preparation & Examinations
- Gynecological cleaning
- Laboratory cell research
- For cell collection in veterinary applications

Ref #	Packaging
816 <i>Non-sterile</i>	10/50 (10 boxes per case, 50 per box)
816 100 <i>Non-sterile</i>	10/100 (10 boxes per case, 100 per box)
816 BULK <i>Non-sterile</i>	500 bulk (10 Bdls. per case, 50 per Bndl)

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com