




Polyester Tipped Applicator

Plastic handle, inherent score mark



Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It is non-absorbent and contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a tapered plastic handle by an aqueous based adhesive. Inherent score mark on handle. Available sterile or non-sterile.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.078" Dia. (1.981 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.085" Dia. (2.159 mm) tapered to 0.078" (1.981 mm)
- 5.91" Length (150.114 mm)

Qualities

- Medical grade quality
- Non-absorbent tip
- Synthetic polyester contains no fatty acids or latent sugars
- Small tip for patient comfort
- Molded handle is semi-flexible
- Made in USA

Applications

- Diagnostic specimen collection
- Recovery of micro-organisms and fluid
- Cell collection for DNA testing
- For use in a wide range of diagnostic tests and laboratory uses

Ref #	Packaging
1000-PD Non-Sterile	10/100 (10 bags per case, 100 swabs per bag)
25-1000 1PD Sterile	10/100/1 (10 boxes of 100 individually wrapped packages, 1 swab p/pkg)



Puritan[®]